

McLaren Print System Order

Order No: 56993
 Order Date: 2020-09-22
 User: Victoria Tijerina
 Phone: 5173031371

Ship Location: Grand Ledge Health Center
 1035 Charlevoix Dr Ste 200
 Grand Ledge , MI 48837

Brochures
 Quantity: 25
 Paragon Dept No: 51015
 Dept Name: McLaren Grand Ledge
 Company Number: 810


Order Total Price: 20.00

Item Number: M-10370
 Item Description: Breastmilk Storage Magnet
 Revision Date: 2014
 Print:
 Paper:
 Size:
 Fold:
 Finish:
 Drill:
 Misc Info: 4x7 Magnet with rounded corners

Breast Milk Storage Guidelines for Healthy Babies

Freshly expressed breast milk	Room Temperature 66-72°F	Refrigerator 32° F	Freezer 4° F <small>(Don't store on door. Keep frozen solid.)</small>
	4 hours	5-7 days	3-4 months 6-12 months or longer in deep freezer (0° F or less)

- Wash hands before pumping.
- Use "oldest" milk first. Label milk with date.
- Thaw milk in refrigerator overnight or hold bottle under warm running water.
- To bring milk up to body temperature, place sealed container in a bowl of warm water for 20 minutes.
- Thawed milk is safe in refrigerator for 24 hours.
- NEVER refreeze milk.
- NEVER microwave breast milk.

Spec Info:  **McLaren** *If you have questions or concerns, call McLaren's Lactation Consultant at (810) 342-4896.*
 FAMILY BIRTHPLACE
 M-10370 (2014)