

McLaren Print System Order

Order No: 83162
 Order Date: 2024-02-14
 User: charnell mcqueen
 Phone: (810) 342-2279

Ship Location: Charnell McQueen
 401 S. Ballenger Hwy
 Flint, MI 48532

Brochures
 Quantity: 100
 Paragon Dept No: 23070
 Dept Name: 7S LDRP
 Company Number: 60

Order Total Price: 80.00

Item Number: M-10370
 Item Description: Breastmilk Storage Magnet
 Revision Date: 2014
 Print:
 Paper:
 Size:
 Fold:
 Finish:
 Drill:
 Poster:
 Misc Info: 4x7 Magnet with rounded corners

Breast Milk Storage Guidelines for Healthy Babies

Freshly expressed breast milk	Room Temperature 66-72°F	Refrigerator 32° F	Freezer 4° F <small>(Don't store on door. Keep frozen solid.)</small>
	4 hours	5-7 days	3-4 months 6-12 months or longer in deep freezer (0° F or less)

- Wash hands before pumping.
- Use "oldest" milk first. Label milk with date.
- Thaw milk in refrigerator overnight or hold bottle under warm running water.
- To bring milk up to body temperature, place sealed container in a bowl of warm water for 20 minutes.
- Thawed milk is safe in refrigerator for 24 hours.
- NEVER refreeze milk.
- NEVER microwave breast milk.

Spec Info:



*If you have questions or concerns, call
 McLaren's Lactation Consultant at (810) 342-4896.*

